

# Two-phase Wittig reaction in a continuous flow microreactor

## Background

A convenient and selective way of forming a carbon-carbon double bond is through the Wittig reaction and its modifications. Although easily performed in batch, the Wittig reaction could be incorporated in a total synthesis using continuous flow. The latter has the added advantage of handling all toxic and corrosive reagents inside a closed system.




Figure 1: Two-phase Wittig reaction scheme

In conventional Wittig reactions, the phosphonium ylide reagent is prepared in batch by treatment of the Wittig salt with a base. However, when using substrates which can withstand alkaline conditions, the ylide can be prepared in situ. In this so-called two-phase Wittig reaction, an aldehyde (1) and a Wittig salt (2) are introduced into the microreactor, where they react to form the corresponding unsaturated Wittig product (3) by treatment with aqueous sodium hydroxide. The Wittig reaction is an example of a very selective reaction, producing no side products. However, the formed double bond can take two configurations: the *cis* (Z) and the *trans* (E) form. The latter is shown in Figure 1.

## Setup and method

### Material

- FlowStart B-200
- B-230 Pump Module
- B-242 Inlet Module
- Basic Quench Microreactor (internal volume  $V_{\mu\text{R}} = 92 \mu\text{L}$ )




Figure 2: FlowStart setup for the two-phase Wittig reaction

### Chemicals

Recommended grade: *pro analysi* (p.a.) or *reagent grade*.

- tert-Butoxycarbonylmethyl phosphonium bromide (Wittig salt)
- benzaldehyde
- sodium hydroxide
- acetic acid
- dichloromethane (DCM)
- anisole (used as internal standard)

### Stock solutions

- A.** Benzaldehyde (102  $\mu\text{L}$ , 1.00 mmol), anisole (109  $\mu\text{L}$ , 1.00 mmol) and the Wittig salt (457 mg, 1.00 mmol) dissolved to a total volume of 10 mL with dichloromethane (corresponding to 0.1 M)
- B.** Sodium hydroxide 0.2 M in water (corresponding to 0.2 M)
- Q.** Acetic acid (1.73 mL, 30.0 mmol) dissolved to a total volume of 10 mL with dichloromethane (corresponding to 3.0 M)

Stock solutions are to be prepared at the beginning of the experiments. Make sure to close the flasks which are used to store the solutions, as some of the components are rather volatile.

### Analysis

Analysis of the reaction mixture is done using gas chromatography. Calibration of the product against the internal standard (anisole) is done using general methodology. For a quick calibration, make 4 samples with a fixed anisole concentration and a varying compound concentration. Analyse these samples and 1) setup a calibration curve of peak area ratio against concentration ratio and/or 2) determine the relative response factor.

### Basic experiment

To get acquainted with the reaction and with flow chemistry in general, a so-called *basic experiment* is performed. This experiment is the Wittig reaction at fixed parameters – a reaction time ( $t_R$ ) of 5.0 min, a temperature of 20°C and a sodium hydroxide molar excess ratio ( $ME_{B/A}$ ) of 2.0. The Q/B molar excess ratio ( $ME_{Q/B}$ ) is set to a fixed value of 20. The target volume of solution A to be collected is 20  $\mu\text{L}$ , all samples are collected in a GC vial containing 500  $\mu\text{L}$  dichloromethane. The used setup can be seen in Figure 2.

The corresponding flow rates can be calculated according to the known equations. After preparation of this experiment, the instructor should check if the calculated flow rates and collection time are correct.

### Procedure

- Prepare solutions A, B and Q
- Fill the three syringes with solutions A, B and Q
- Slide the microreactor into the holder and connect inlet and outlet tubing
- Connect the inlet tubing to the corresponding syringes, and place the syringes on the pumps
- Set the right flow rates and press start
- Stabilise for 15 minutes
- Collect your sample for the calculated time
- Analyse your sample using gas chromatography and calculate substrate conversion from the calibration curve or relative response factor
- Rinse the *FlowStart* system by purging the tubing and microreactor with acetone

- Empty, clean and dry the syringes afterwards

**Note:** Make sure to close the vial after collecting. This is done because some of the reaction components are rather volatile and readily evaporate from the vial.

## Questions

1. **Preparation of the experiment:**
  - a. Roughly calculate the cost of the experiment from the prices of the chemicals. In other words, calculate the price (e.g. per gram) of the product.
  - b. Find the safety aspects (including R/S values) of the used chemicals.
2. **Q:** Find the reaction mechanism for the Wittig reaction, and show the essential (sequential) steps from benzaldehyde and phosphonium ylide to Wittig product.
3. **Q:** What advantages in performing the Wittig reaction in continuous flow can you think of? Also, can you think of any disadvantages?
4. **Q:** Because the sodium hydroxide in water solution does not mix with the dichloromethane solution, this Wittig reaction is carried out in a two-phase system (comparable to an 'oil-in-water' system). You can observe this in the microreactor by the formation of so-called *plugs*. The actual reaction takes place at the *interface* (area of contact) of these plugs. Can you think of advantages to this plug-flow reaction? Also, what are the disadvantages? What possible solutions could circumvent this two-phase system, and what are its (dis)advantages?

